VR11T

PK CMV-PA SYSTEM Passive Particle Agglutination Test for Detection

I. INTENDED USE

particle agglutination assay intended for the recipients. qualitative detection of IgG and IgM antibodies to cytomegalovirus (CMV) in human EDTA A variety of methods have been developed to detect antibodies to CMV including indirect plasma and serum from blood donors using the Beckman Coulter PK7300 and/or PK7400 Automated Microplate Systems. This test is not intended for diagnostic use.

II. SUMMARY OF TEST

Cytomegalovirus (CMV) is a double-stranded DNA virus with physicochemical characteristics variety of protein antinens to the surfaces of Serologic surveys have shown that CMV infection is worldwide in distribution, with indefinitely following initial infection, or it may emerge from time to time to cause an active infection.⁵ Pregnant women can transmit the CMV virus to the fetus, resulting in congenital liver, spleen, and/or CNS disease in the unborn Automation has enhanced the value of the immunosuppressed patients such as organ or PK7400 Automated Microplate Systems. transplant recipients and in patients harboring human immunodeficiency virus (HIV).^{8,9} III. PRINCIPLE OF PROCEDURE The PK CMV-PA System uses gelatin particles

The transfusion or transplant of CMV- coated with cytomegalovirus antigens to detect seropositive blood or organs may cause a variety IgG and IgM antibodies to CMV in human recipients. Since CMV infections are frequently material is diluted with SAMPLE DILUENT

Catalog No. PA2300 November, 2023

of Total Cytomegalovirus Antibodies using Beckman Coulter PK Systems

toward reducing cytomegalovirus infection in The PK CMV-PA System is a passive immunocompromised transfusion and transplant

to detect antibodies to CMV including indirect hemagglutination assay (IHA), indirect fluorescent assay (IFA), anticomplement immunofluorescence (ACIF), enzyme immunoassay (EIA), or passive latex agglutination (PLA).

common to members of the herpesvirus family.1 tannic acid-treated sheep erythrocytes and demonstrated hemagglutination in the presence antibody prevalence in adults in the range of these methods are still widely used today despite of the corresponding antibodies. Variations of 20-82%, ^{2-4,18} The majority of CMV infections are subclinical or associated with nonspecific oddrage these problems associated with biological carriers. To are subclinical or associated with nonspecific illness. The virus may remain in a latent state indefinitive following initial inclusion in a latent state

child.⁶ The congenital effects of mother-to- indirect particle agglutination test by significantly fetus CMV transmission may be more severe in reducing the amount of time and labor needed to those cases where the mother has acquired the perform the assay. The PK CMV-PA System has primary infection during early pregnancy than been developed to provide an indirect particle in maternal cases of reactivated disease.⁷ CMV agglutination CMV assay using uniform reagents infection often induces life-threatening conditions which are stable, easy to handle, and suitable for such as pneumonia, fever, and hepatitis among use on the BECKMAN COULTER PK7300 and/

of clinical abnormalities in immunocompromised serum and plasma. The test sample or control transmitted through organ transplants and blood and then mixed with the sensitized particles transfusions,^{2-4,10-12} the screening of blood in a terraced microplate well. During the donors for CMV antibodies is an important step incubation, the particles settle in the terraced

microplate well. Antibody to CMV will bind to 2) The microplates must be clean and in the antigen-sensitized particles during this good condition before use. Damaged incubation. Particles with bound antibody will form addlutination, which are visible as a homogeneous blue layer of gelatin particles. When antibodies to CMV are not present, sensitization and subsequent addlutination does not occur. Particles without bound antibody fall freely to the center of the well and visually appear as a compact dense blue button surrounded by a clear zone.

The PK7300 and/or PK7400 instrument will read the settling patterns of particles in each well based on the threshold settings chosen for the reagent. The PK7300 and/or PK7400 determines the presence or absence of antibodies to CMV using a CCD (charged coupled device) camera, which captures the well image allowing differentiation of agglutinated and unagglutinated 3) If positive control samples repeatedly patterns.

IV.REAGENTS

The PK CMV-PA System is available in a kit sufficient to perform 2500 tests. Store reagents at 2-8°C. DO NOT FREEZE.

SENSITIZED PARTICLES - 10 vials containing gelatin particles colored with blue dye, sensitized with cytomegalovirus antigens and then lyophilized. Each vial must be reconstituted 5) Sodium azide is added to the reagents as a with 6.0 mL RECONSTITUTING SOLUTION. Reconstituted particles contain 0.15% sodium azide. Each vial is sufficient for 250 tests on the PK7300 and on the PK7400.

RECONSTITUTING SOLUTION -1 bottle, 70 mL. Phosphate buffered saline containing 0.10% sodium azide. For the reconstitution of 6) Visible signs of microbial growth or gross SENSITIZED PARTICLES.

SAMPLE DILUENT - 3 bottles, 300 mL each. Proprietary solution containing phosphate buffered saline, normal rabbit serum and 0.10% 7) Handle all specimens, human-based sodium azide.

V. WARNING AND PRECAUTIONS

PK CMV-PA System is for in vitro diagnostic use. 1) Avoid contamination of reagents or specimens with saliva which can cause 8) Do not eat, drink or smoke in areas where indistinguishable agglutination patterns. Do specimens, human-based reagents and not mouth pipette any reagents.

Sample ID

PK7300

NR

Expected Result

TABLE 6. REFERENCE METHOD (PK7300) COMPARISON TO PK7400 - PLASMA: ALL LOTS/SITES

RNRTotalAgreementTotalAgreementConfidence
bound148501485PPA14851485100.00%99.75%

2 1885 1887 **NPA** 1885 1887 99.89% 99.62%

 # correct Result
 % Correct

 Channel 11
 Channel 12
 Result

100%

100%

100%

100%

100%

100%

100%

100%

100%

100%

Total 1487 1885 3372 OPA 3370 3372 99.94% 99.79%

ABLE 8. REPRODUCIBILITY OF THE PK CMV-PA SYSTEM ON THE PK7400 FOR ALL LOTS/SITES

 2400016401
 Pos
 30/30
 30/30

 2400016402
 Pos
 30/30
 30/30

 2600016403
 Pos
 30/30
 30/30

 2600016404
 Pos
 30/30
 30/30

 2600016404
 Pos
 30/30
 30/30

 9245209
 Neg
 30/30
 30/30

 2400016406
 Pos
 30/30
 30/30

103459900 Neg 30/30 30/30

 1034550100
 Neg
 30/30
 30/30

 1034550400
 Neg
 30/30
 30/30

 1034550500
 Neg
 30/30
 30/30

 1034551000
 Neg
 30/30
 30/30

 1034549600
 Pos
 30/30
 30/30

plate terraces or improper washing of the microplates resulting in protein buildup or debris in the terraces can adversely affect test results. For instance, the homogenous laver of agglutinated particles in a positive reaction can be disrupted and fold over onto itself. An analogy would be the folding over of the sides of an omelet. This phenomenon can result from excessive vibration, protein buildup in the terraces of the microplate or physical damage to the microplate terraces, and is readily apparent during plate review. The recommended microplate maintenance procedures can be found in the Beckman Coulter PK7300 User's Guide and the PK7400 Instructions for Use.

- test negative, excessive instrument vibration is a potential cause. When control material repeatedly fails to perform as expected, contact Beckman Coulter Immunohematology Technical Services at 800-447-5852.
- 4) Avoid freezing of PK CMV-PA reagents and reconstituted SENSITIZED PARTICLES.
- bacteriostatic agent. Sodium azide has been reported to form explosive lead and copper azides in laboratory plumbing. To prevent azide build-up, flush with large volumes of water if solutions containing azide are disposed of in the sink.
- turbidity in the reagent may indicate degradation and warrant discontinuance of
- reagents and controls as if potentially infectious. Refer to the Center for Disease Control guidelines for handling biological materials.1
- controls are handled.

- donor samples.
- most cases.

- at 2-8°C.
- reconstituted expiry should be recorded on the reagent vials. 4) After the reconstitution period, gently swirl (DO NOT VORTEX) the reagent to assure thorough resuspension. Place the

All Sites/Lots	Tri	al PK74	400	Statistical Summary		Statistical Summary				
Reference PK7300	R	NR	Total		Agreement	Total	Rate of Agreement (%)	Lower 95% confidence bound		
R	1690	2	1692	PPA	1690	1692	99.88%	99.57%		
NR	6	1991	1997	NPA	1991	1997	99.70%	99.35%		
Total	1696	1993	3689	OPA	3681	3689	99.78%	99.57%		
TABLE 7. REFER		NETHO	D (PK73	00) CO	MPARISON TO) PK740	0 - SERUM: AL	L LOTS/SITES		
All Sites/Lots	Trial PK7400				Statistical Summary					
Poforonco							Rate of	Lower 95%		

Sample Type	True Result Positive	Incorrect Negative Result	Sensitivity	Lower 95% confidence bound
EDTA Plasma	1690	0	1690/1690 100.0%	99.78%
Serum	1485	1	1485/1486 99.93%	99.63%

Sample EDTA Plas Serun

method comparison testing to the "true" res obtained on the PK7300. Discordant sample were tested by two other methods (DiaSorin summarized in Tables 10 through 13. LIAISON assay and Immucor Capture-CMV).

PK-CMV-PA System				PLASMA				SERUM		
PK7200	PK	7300		Initial		Repeat	In	Initial		Repeat
Neg	N	eg		3340		3365	14	455		1462
Neg	P	os	T	19		2		1		0
Pos	N	eg	T	15		1		11		2
Pos	P	os		4020		4026	1:	516		1519
Tot	als		T	7394		7394	29	983		2983
% CONCO	ORDAN	NCE		99.50%		99.9%	99	.6%		99.9%
TABLE 2. RA	TE OF	AGRE	EEM	IENT FOR	SAI	MPLE AG	E SUBS	ET C		E PK7300
Test Da	у			iber In ement		Rate o Agreemer			Lower fidenc	95% e bound
Initial		1	1143	13/1154		99.05%			98.43%	
Day 4		ſ	1144	4/1154		99.13%			98.53%	
Day 6 11		1146	6/1154 99.31%		%	98.75%		5%		
ABLE 3. RE	PROD	UCIB	ILIT	Y OF TH	E PK	CMV-PA	SYSTE	ЕМ О	N THE	PK7300
				N =	D	AY 1-3	DA	Y 4	1	DAY 6
Known R	eactive)		81		81	7	6		81
Known Non	-React	ive		12		11*	11	*		11*
1 sample te	sted fa	lse po	sitiv	/e on day	1-3 a	and 6.				
TABLE 4.	SENS		ry c	OF THE P	K CI	/IV-PA SY	STEM	он т	HE P	(7300
Sample Typ	e I	Correc Resul ositiv	t	Incorre Negati Resul	ve	Sen	sitivity		(95%	itivity % % lower Bound)
Serum		1517		4		1517/15	21 99.7	4%	99	.40%
EDTA Plasm	a	4026		0		4026/40	26 100.	0% 99.93%		.93%
TABLE 5.	SPEC	IFICIT	ry c	OF THE P	кс	/IV-PA SY	STEM	ON 1	HE P	(7300
Sample Typ	e I	Correc Resul egativ	t	Incorrect Positive Result		Spe	cificity		(95%	ificity % 6 lower Bound)
Serum		1462		0		1462/14	62 100.	0%	99	.80%

EDTA Plasma 3355 13 3355/3368 99.61% 99.39%

TABLE 1. INITIAL AND REPEAT CMV RESULTS WITH PLASMA AND SERUM

sult	Best two out of three results was considered the
ples	"true" result after additional testing. Results are
orin	summarized in Tables 10 through 12

9) Clean pipettes should be used to reconstitute lot number may be pooled following all reagents. Clean glass or plastic containers completion of the reconstitution period. The should be used for pooling reagents from the

10) Positive and negative control materials should be handled in the same fashion as 6) SENSITIZED PARTICLES from one lot

11) Inadequate adherence to the package insert can result in erroneous results. 7) The SAMPLE DILUENT AND

12) Carrvover between samples has been detected in some donor samples with high titers of CMV antibodies. The PK7300 and Note: All reagents should be brought to room PK7400 require the use of Cleaning Solution (no preparation required) in the Cleaning and use Solution Tank, which eliminates carry-over in VII. STORAGE

VI. REAGENT PREPARATION

1) Reconstitute, as needed, each vial of SENSITIZED PARTICLES with 6 mL of RECONSTITUTING SOLUTION. Replace the stopper and invert a few times to assure thorough mixing. Prior to use, allow the 2) Store the PK CMV-PA System at 2-8°C. DO reagent to reconstitute for 3-10 hours at NOT FREEZE. room temperature (15-30°C).

3) The PK CMV-PA System should not be used **Note: The reagent reconstitution time** after the expiration date which is printed on has been updated from a minimum of 30 the outside of the package. minutes to 3-10 hours.

4) Store reconstituted SENSITIZED 2) Reconstituted particles are stable for 7 days PARTICLES at 2-8°C. DO NOT FREEZE. SENSITIZED PARTICLES are stable for 7 days after reconstitution, when stored at 3) The date of reconstitution and the

in the mixture.

another lot number.

expiration date.

technique.

number should not be mixed with those of

RECONSTITUTING SOLUTION are not

matrixed to the SENSITIZED PARTICLES lot.

temperature (15-30°C) before reconstitution

Note: Reagents should not be used after the

1) Unused particle suspension should be

returned to the original vial, using aseptic

5) Visible signs of microbial growth or gross turbidity in the reagents may indicate degradation and warrant discontinuance of

reconstituted well mixed sensitized particles into a clean dry PK reagent vial. Place the PREPARATION

PK vial into the PK reagent tray and place Plasma (EDTA) and serum samples, obtained Plasma (EDTA) and serum samples, obtained through standard collection procedures are suitable for this assay. The performance of this assay has not been established with plasma samples employing heparin as the anticoagulant, serum samples collected with serum separator 5) SENSITIZED PARTICLES from the same tubes, heat-treated samples, or neonatal

mixture is stable for 7 days from the earliest samples, pleural fluid, saliva, or nonhuman SYSTEM: reconstitution date of the particles contained samples.

Prior to analysis on the PK7300 and/or PK7400, samples should be adequately centrifuged to ensure that the plasma or serum is free from MATERIALS REQUIRED BUT NOT PROVIDED: particulate matter. If erythrocytes or other visible -Beckman Coulter microplates with a 5 µm well components are contained in the sample, terraces remove by centrifugation to prevent interference -Pipetting device capable of delivering 6.0 ml with the test results. The PK7300 User's Guide -BECKMAN COULTER PK7300 and/or and/or the PK7400 Instructions for Use requires BECKMAN COULTER PK7400 centrifugation of samples within 10 hours of -PK CMV-PA System Controls centrifugation of samples within 10 hours of analysis and centrifugation for a minimum of 10 minutes at 1000 x g. These requirements exist for the purpose of optimizing red cell sampling. Therefore, plasma or serum samples tested do not need to comply with these requirements as long as the plasma or serum is free from particulate material. Samples exhibiting cross linemia, hemolysis or icterus may be gross lipemia, hemolysis or icterus may be the PK7300 User's Guide or Chapter 3 of the compromised and may require alternative PK7400 Instructions for Use. testing.

Test parameters and recommended thresholds for the PK7300 and the PK7400 have been established based upon application development with characterized samples. Working files for the PK CMV-PA test are shown below for the PK7300 and PK7400. Good laboratory practice dictates that each laboratory validate the operating parameters. may be stored frozen at<-20°C if testing is to All reagents, diluents, and specimens should be exceed 14 days after collection. Samples should at room temperature (15-30°C) prior to analysis. be well mixed after thawing. Repeated freeze/ thaw cycles should be avoided. Improper storage of specimens may result in variable settling patterns yielding false positive or indeterminate

When shipping specimens, they should be packaged in compliance with applicable federal, state and local regulations covering the transport of clinical specimens and etiologic agents. Specimens may be shipped at either ambient. refrigerated (2-8°C) on wet ice, or frozen (-10°C or colder) on dry ice.

samples. In addition, the performance of this IX. MATERIALS assay has not been established with cadaveric MATERIALS PROVIDED IN THE PK CMV-PA -RECONSTITUTING SOLUTION

-SENSITIZED PARTICLES

	Threshold Settings for the PK7300 and PK7400								
	F	P/C	S	РС	L	A			
	(+) Limit	(-) Limit	Low	High	(+) Limit	(-) Limit	LIA Selection	BG/C	
PK7300	35	20	14	14	200	90	5	Middle	
PK7400	33	19	18	18	200	90	5	Middle	
Dynamic Range Settings for the PK7300 and PK7400 P C LIA SPC									
		High	Low	Higl	n Low	High	Low	High	
	Low					450	N1/A	N/A	
PK7300	50	99	10	99	0	450	N/A	IN/A	

PK7300 PARAMETERS VOLUME/SETTING

ample/Diluent Ratio	8.9	Sample/Dilu
iluted Sample Volume	25 <i>µ</i> L	Diluted San
eagent Volume	23 µL	Reagent Vo
eagent Name	CMV	Reagent Na
hannel Designation	(11-12)	Channel De
ecision Logic	Positive/ Negative	Decision Lo
emperature Setting	28°C	Temperatur
ncubation Time	60 minutes	Incubation
late Well	5 <i>µ</i> m	Plate Well

NOTE: CHANGES MADE TO THE PK7300 and/or PK7400 ARE AUTOMATICALLY SAVED TO THE HARD DRIVE WHEN THEY ARE MADE. IT IS RECOMMENDED THAT WHEN ANY CHANGES ARE MADE, THEY ALSO BE SAVED TO AN EXTERNAL STORAGE MEDIA.

XI. QUALITY CONTROL

NONREACTIVE CONTROLS should be tested making sure that the volume of the controls is at the beginning and end of each batch of sufficient for adequate instrument sampling (> samples assayed, after the addition of reagents 1.5 mL). When control material repeatedly fails to

package insert for complete details regarding 447-5852. this material. Additional quality control testing XII. INTERPRETATION may be performed by the user including other The PK7300 and PK7400 will interpret the well- characterized specimens or referenced settling patterns of particles in each well sera.

User's Guide or the PK7400 Instructions for Use the PK7400 Instructions for Use, for complete using the reactive and nonreactive controls as details of the analyzer's interpretation of specimens. The reactive control should produce reactions. a positive reaction and the negative control The PK7300 and PK7400 stores the reaction should produce a negative reaction with the test. If appropriate results are not obtained with

TABLE 9. REPEATABILITY OF THE PK CMV-PA SYSTEM ON THE PK7400 FOR ALL LOTS/SITES

Commis ID	Expected	# Correc	% Correct	
Sample ID	Result	Channel 11	Channel 12	Result
2400016401	Pos	15/15	15/15	100%
2400016402	Pos	15/15	15/15	100%
2600016403	Pos	15/15	15/15	100%
2600016404	Pos	15/15	15/15	100%
9245209	Neg	15/15	15/15	100%
2400016406	Pos	15/15	15/15	100%
103459900	Neg	15/15	15/15	100%
1034550100	Neg	15/15	15/15	100%
1034550400	Neg	15/15	15/15	100%
1034550500	Neg	15/15	15/15	100%
1034551000	Neg	15/15	15/15	100%
1034549600	Pos	15/15	15/15	100%

TABLE 10. SENSITIVITY OF THE PK CMV-PA SYSTEM ON THE PK7300

TABLE 11 SPECIFICITY OF THE PK CMV-PA SYSTEM ON THE PK7300

е Туре	True Result Negative	Incorrect Positive Result	Specificity	Lower 95% confidence bound	
Plasma	1997	2	1997/1999 99.90%	99.64%	
um	1886	0	1886/1886 100.0%	99.80%	

TABLE 12. SENSITIVITY OF THE PK CMV-PA SYSTEM ON THE PK7400

Sample Type	True Result Positive	Incorrect Negative Result	Sensitivity	Lower 95% confidence bound
EDTA Plasma	1690	0	1690/1690 100.0%	99.78%
Serum	1486	0	1486/1486 100.0%	99.75%

TABLE 13. SPECIFICITY OF THE PK CMV-PA SYSTEM ON THE PK7400

Sample Type	True Result Negative	Incorrect Positive Result	Specificity	Lower 95% confidence bound				
EDTA Plasma	1993	6	1993/1999 99.70%	99.35%				
Serum	1885	1	1885/1886 99.95%	99.70%				

XVII. REFERENCES

- 1. Griffiths PD, Grundy JE. Molecular biology and immunology of cytomegalovirus. Biochem. J. 1987; 241:313-324.
- 2. Perham TGM, Caul EO, Conway PJ, Mott MG. Cytomegalovirus infection in 8. Ho M. Viral infections after transplantation blood donors - a prospective study. Brit. J. in man. Arch. Virol. 1977; 55:1-24. Haem. 1971; 20:307-320.
- 3. Prince AM, Szmuness W, Millian following renal transplantation. Rev. Inf. SJ, David DS. A serologic study of Dis. 1981; 3:1151. cytomegalovirus infections associated with 10. Beebe JL. Cytomegalovirus: disease and blood transfusions. N. Engl. J. Med. 1971; 284:1125-1131.
- 4. Pass MA et al. Evaluation of a walking 11. Yeager AS et al. Prevention of transfusiondonor blood transfusion program in an intensive care nursery. J. Pedial. 1976;89: 646-651.
- 5. Drew WL. Cytomegalovirus: an overview. 12. Adler SP. Transfusion associated Syva. Monitor. 1988;6(2):1-5.
- 6. Yeager AS. Transfusion-acquired Dis. 1983; 5:977-993. cytomegalovirus infections in newborn infants. Am. J. Dis. Child. 1974;128: 478-13. Boyden SV. The adsorption of protein on erythrocytes treated with tannic acid

- Stagno S et al. Congenita cytomegalovirus infection: the relative importance of primary and recurrent maternal infection. N. Engl. J. Med. 1982; 306:945-949.
- 9. Glenn J. Cytomegalovirus infections
- diagnosis. Lab Management 1978; 16:15-

acquired cytomegalovirus infections in newborn infants. J. Pediat. 1981; 98:281-287.

cytomegalovirus infections. Rev. Infect.

14. Kobavashi N et al. Establishment of the technology of a new specific diagnostic test for syphilis using Treponema antigen (MDA-TP). STD 1985; 66:61-73.

193:107-120.

and subsequent hemagglutination by

antiprotein sera. J. Exp. Med. 1951;

- 15. Deguchi M et al. Measurement of antiponema pallidum antibodies using a new artificial carrier, high-density composit particles (HDP). Med and Pharm 1991; 26:333-339.
- 16. Yoshida T et al. Evaluation of passive particle agglutination test for antibody to human immunodeficiency virus. J. Clin Microbiol. 1987: 1433-1437.
- 17. US Department of Health and Human Services. Biosafety in microbiological laboratories. HHS Publication (NH) 1-1112, Washington US Government Printing office, December 2009.
- Killian CS. Prevalence of cytomegalovirus antibodies in blood donors. Lab medicine 1989; Be:103-105.

VOLUME/SETTING

ent Ratio	8.9
ple Volume	25 <i>µ</i> L
ume	23 µL
me	CMV
signation	(11-12)
gic	Positive/ Negative
Setting	30°C
ime	60 minutes
	5 <i>µ</i> m

the controls, all assay results within that batch The PK CMV-PA System REACTIVE and are invalid and must be retested. Repeat testing and after interruption or delays in processing. perform as expected, contact Beckman Coulter Refer to the PK CMV-PA System Controls Immunohematology Technical Services at 800-

based on the threshold settings chosen in the Perform testing as described in the PK7300 parameter file. See the PK7300 User's Guide or

of the sides of an omelet. This phenomenon can for the reaction. A negative SPC value together in the terraces of the microplate or physical cause the channel result to be indeterminate. Anage to the microplate terraces, and is readily apparent during plate review. The recommended microplate maintenance procedures can be found in the PK7200 Hear's Guide or the found in the PK7300 User's Guide or the to CMV, indicating that the individual has not PK7400 Instructions for Use.

A complete description of plate inspection and results review is contained in Section C of the PK7300 User's Guide or Chapter 2 of the CMV by the criteria of the PK CMV-PA System. must be performed on any samples for which visual and analyzer interpretations do not agree. Refer to Section C of the PK7300 User's Guide or Chapter 2 of the PK7400 Instructions for Use.

A sample reported as indeterminate (?) on The presence or absence of antibody to initial screening may be considered reactive cytomegalovirus is determined by the PK7300 by the criteria of the PK CMV-PA System, may and/or PK7400 using a CCD camera which be repeated in duplicate on the PK analyzer or analyzes the well image and can differentiate tested by an alternative method. agglutinated and unagglutinated patterns. If an initially indeterminate sample is repeated agglutinated and unagglutinated personal The PK7300 and/or PK7400 employs three assessment parameters for each microplate well containing PK CMV-PA System reagent and test

- · SPC Sharpness of the edge of the button reactive for antibodies to CMV by the criteria LIA Quantity of particles in the center of the OK CMV-PA System. If upon repeat
- P/C Ratio of the average light transmittance of the peripheral and central values
- of the peripheral and central values CMV- PA System. The parameters SPC, LIA and P/ C are compared to programmable thresholds to obtain an interpretation (+, 2) for each reaction only those samples which test negative on initial screening or in both duplicate retests should be

an interpretation (+,-,?) for each reaction.

considered negative for antibodies to CMV for

purposes of transfusion.

be performed either manually (within 30 minutes) The most important parameter resulting from or on-line (no time limit). Visually, a reactive test the image analysis system is SPC. If the SPC is a homogenous layer of particles. A nonreactive is determined positive, then either a positive test would result in a compact, dense button or indeterminate LIA or P/C value will result in surrounded by a clear zone. Plate review should an overall positive result interpretation for the include inspection of the reactions for abnormal reaction. A positive SPC value together with settling patterns or for any sample for which a negative value for either the LIA or P/C will visual and analyzer interpretations do not agree. cause the channel result to be indeterminate. If Under certain circumstances, this homogenous the SPC is determined negative, then either a layer of particles can be disrupted and fold over negative or indeterminate LIA or P/C value will onto itself. An analogy would be the folding over result in an overall negative result interpretation

been infected with cytomegalovirus.

screening is considered reactive for antibodies to

The presence of antibodies indicates previous

or current infection. Individuals with antibodies to

the specimen is to be interpreted as repeatedly

testing both duplicate results are nonreactive,

result from excessive vibration, protein buildup with a positive value for either the LIA or P/C will

between chronic and acute CMV infections. This product is only for use in screening blood donors and has not been evaluated as a Sensitivity and Specificity for the PK CMVdiagnostic test for CMV outside the blood bank PA System when tested on the PK7300 was

XV. EXPECTED RESULTS

XVI. SPECIFIC PERFORMANCE

CHARACTERISTICS

summarized in Table 1.

PK7300

XIV. LIMITATIONS OF THE PROCEDURE

circulating antibodies to cytomegalovirus. It has

been shown to be safe and effective for the large block denors when used in the PK CMV-PA System was tested with a scale screening of blood donors when used in accordance with instructions provided. Donors in PK7300 from individuals demonstrating reactivity the earliest stages of infection may not contain detectable levels of CMV antibody. The PK CMV-PA System is not intended to distinguish between chronic and acute CMV infections cross reactivity.

The PK CMV-PA System is used to detect are summarized in Tables 2 and 3.

Several studies have shown the expected field trial testing to the "true" result obtained on

incidence of CMV antibodies in various the PK7200. Discordant samples were tested

populations. In a recent study of 250 random by two other methods (EIA and Capture-CMV).

blood donors, 50% were positive for CMV Best two out of three results was considered the

performance characteristics of the PK CMV-PA The performance of the PK CMV-PA System

System in blood donors were evaluated in two was evaluated on the PK7400 by comparing

sites on the PK7300. Evaluation of 2020 blood to the PK7300 reference results. Testing was

donor samples demonstrate CMV seropositivity performed at three geographically distinct blood

of 43.3%, 50.7%, for Sites 1 and 2, respectively. centers. A total of 3372 serum samples and 3689

The performance of the PK CMV-PA System on the PK7400 was evaluated at three sites

was evaluated on the PK7300 by comparing using three lots (one lot per site) by testing

to the PK7200 reference results. Testing was 12 plasma and serum samples. Out of the 12

performed at two geographically distinct blood samples, 6 were known reactive and 6 known

centers and also at a Beckman Coulter facility. A non-reactive samples. All samples were tested

total of 2983 serum samples and 7394 plasma in duplicate, in two different runs per testing day

(EDTA) samples were tested. Results are (at least 2 hours apart), over a minimum of 5

the PK7300 was evaluated by testing a subset, Sensitivity and Specificity for the PK CMV-

1154 plasma samples, of the total samples tested PA System when tested on the PK7400 was

in the study. Included in the subset were 27 known determined by comparing the 3372 serum

reactive and 4 known non-reactive samples. All samples and the 3689 plasma samples from

 Reproducibility
 Tables 6 and 9.

 The reproducibility of the PK CMV-PA System on
 Sensitivity and Specificity

Reproducibility

Tables 8 and 9.

showing CMV antibody prevalence ranging summarized in Tables 4 and 5.

from 20-82%.^{2-4,18} Expected values may vary with age, sex and geographic location. The

antibodies.³ This study supports earlier studies "true" result after additional testing. Results are

determined by comparing the 2983 serum

samples and the 7394 plasma samples from

plasma (EDTA) samples were tested. Results

The reproducibility of the PK CMV-PA System

nonconsecutive days. Results are summarized in

are summarized in Tables 6 and 7.

samples were tested on days 1-3, 4 and 6. Results

MANUFACTURED BY: Fuiirebio Inc. 1-8-1 Akasaka, Minato-ku Tokvo 107-0052, Japan +81-3-6279-0899 DISTRIBUTED BY: BECKMAN COULTER, INC 250 S. Kraemer Blvd. Brea. CA 92821 Customer Service: 800-526-3821 Technical Support: 800-447-5852 Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries. Made in China

